



## AQ Medical 350 Vaccine Refrigerator

The AQ Medical 350 vaccine refrigerator is purpose-built for vaccine storage and follows the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines.

It was designed with low noise output and consistent temperature readings in mind, and maintains a stable temperature between the recommended +2°C to +8°C. The unit holds safety features such as visual and audible alarms in the event that temperatures fluctuate outside the set range.

All AQ Medical refrigerators are accredited by the Quality Care Pharmacy Program (QCPP), a quality assurance program that was developed by the Pharmacy Guild of Australia in consultation with the Pharmaceutical Society of Australia.

AQ Medical units are also backed by our extensive Avem Quirks service network. With qualified technicians and warehousing facilities across Australia, we can deliver responsive service and exceptional support on a national level.

The AQ Medical 350 vaccine refrigerator is registered with GEMS and has a 7 Star Rating, one of the most energy efficient vaccine refrigerators on the Australian Market. Registration Number: ARC1554



Distributed by:  
AMA Medical Products  
[amamedicalproducts.com.au](http://amamedicalproducts.com.au)  
[sales@amamedicalproducts.com.au](mailto:sales@amamedicalproducts.com.au)  
1800 626 292



# AQ Medical 350 Vaccine Refrigerator

## Key Features

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8° C in line with national 'Strive for 5' storage guidelines
- Audible and visual auto resetting alarms
- Triple glazed glass
- Adjustable shelves
- Lockable door
- Low energy consumption
- Low noise output
- Dedicated technical service line

## Optional Extras

- Uninterruptible Power Supply unit with back up battery
- AQ Box cloud temperature monitoring system
- Comprehensive programmed maintenance

## Product Specifications

|                             |                |
|-----------------------------|----------------|
| Dimensions                  | Height 1881 mm |
|                             | Width 580 mm   |
|                             | Depth 584 mm   |
| Gross Weight                | 89 kg          |
| Gross Capacity              | 350 L          |
| Temperature Range           | +2°C to +8°C   |
| Number of Shelves           | 6              |
| Power Consumption/24 hrs.   | 3.3 kWh        |
| Refrigerant                 | R600a          |
| Warranty (parts and labour) | 24 months      |



Distributed by:  
AMA Medical Products  
[amamedicalproducts.com.au](http://amamedicalproducts.com.au)  
[sales@amamedicalproducts.com.au](mailto:sales@amamedicalproducts.com.au)  
1800 626 292

